



Anal irrigation

Peristeen



Peristeen Anal Irrigation

– makes it easier to empty your bowels effectively and securely.

Anal irrigation is a recognised and effective method of irrigating your bowels. If used on a regular basis, **Peristeen** Anal Irrigation can help prevent constipation and faecal incontinence.

A complete system

Peristeen Anal Irrigation contains everything you need to irrigate your bowels. It consists of a control unit with a pump, a water bag and a rectal catheter.

Simple handling

Many details on the **Peristeen** Anal Irrigation have been specially developed for people with low dexterity.

Ergonomic design

The control unit is easy to handle and its self-explanatory symbols guide the user through the procedure.

The pump is easy to use either with the hand, the underarm or the elbow.

The pump may be fastened to the thigh with the strap provided.

The water bag is designed to stand on the floor and can easily be moved around by using the unique large handle.

The rectal catheter is smooth, small and discrete which makes insertion into the rectum very easy.

The soft balloon ensures that the rectal catheter is fixed inside the bowel so both hands are free during the irrigation.



Peristeen Anal Irrigation	
Product code	Description
29121	System
29122	Accessory unit
29123	Rectal catheter
29124	Strap
29125	Tube with blue connectors



Coloplast A/S
Holtedam 1
DK-3050 Humlebæk
Denmark

www.coloplast.com