

# Security is...

The **SpeediCath** product assortment

## SpeediCath

Hydrophilic coating, ready-to-use, PVC free

CH	Male 40 cm	Tiemann Male 40 cm	Female 20 cm	Paediatric 20 cm	Boy 30 cm
• 6			28506	28706	28606
• 8	28408		28508	28708	28608
• 10	28410	28490	28510	28710	28610
• 12	28412	28492	28512		28612
• 14	28414	28494	28514		
• 16	28416		28516		
• 18	28418				

## SpeediCath Complete

Easy and hygienic catheterisation anytime and anywhere

CH	Male 40 cm	Female 20 cm
• 10	28460	28430
• 12	28462	28432
• 14	28464	28434

## SpeediCath Compact

Discreet and tailor-made for women

CH	Female 7 cm	Girl 7 cm
• 8	28578	28578
• 10	28580	
• 12	28582	
• 14	28584	



...freedom to be spontaneous



Coloplast and SpeediCath are registered trademarks of Coloplast A/S. © 2007.09. All rights reserved Coloplast A/S, 3050 Humlebeak, Denmark.



Coloplast A/S  
Holtedam 1  
DK-3050 Humlebeak  
Denmark

[www.coloplast.com](http://www.coloplast.com)

# SpeediCath



# Security is freedom

## A small revolution for women who self-catheterise

Many women who self-catheterise have long wished for a more convenient, discreet catheter; a catheter that can be used easily and safely both inside and outside the home environment. **SpeediCath Compact** is Coloplast's solution to this need.

**SpeediCath Compact** is a hydrophilic-coated, sterile intermittent catheter for single use. It is packaged in saline solution and is ready-to-use right out of the pack.

## Tailor-made for women

The female urethra on average is just 3-5 cm long – significantly shorter than the male urethra. With a catheter length of 7 cm, **SpeediCath Compact** is tailor-made for the typical female anatomy. In a clinical study, **SpeediCath Compact** was documented to be at least as efficient at emptying the bladder<sup>1)</sup> as a conventional female catheter with a length of 15-20 cm.

## Ready-to-use

**SpeediCath Compact** eliminates concerns over catheter preparation and quality of the hydration source. The catheter comes pre-hydrated in a ready-to-use package, it could not be easier to use. Simply break the seal and pull open the telescopic packaging - the catheter is ready to insert.

## Easy and hygienic to insert

**SpeediCath Compact** supports hygienic non-touch insertion. The catheter is fitted with an easy-grip handle that keeps fingers safely away from the sterile part of the catheter and improves user control during insertion thereby reducing catheter contact with the septic environment.

## Discreet

**SpeediCath Compact** provides added discretion for the catheter user. When un-opened, the catheter package is no larger than a lipstick. This ingenious packaging design makes the catheter extremely discreet – without compromising performance.<sup>1)</sup>



## Free of PVC and phthalates

**SpeediCath Compact** is 100% PVC-free and therefore does not contain DEHP or other phthalates. The risk of harmful consequences to the environment is thereby reduced compared to PVC-containing catheters.

<sup>1</sup> Biering-Sørensen et al. Residual urine after intermittent catheterization in females using two different catheters. Scand J Urol Nephrol. 2007;41(4):341-5.