


Biatain® Silicone
Soft, effective and easy to use



Biatain® Silicone

The superior absorbency of foam and the gentle touch of silicone




Easy application for safe and correct use

Biatain Silicone has a unique 3-part opening for non-touch application.

This allows you to open and apply the dressing to a wound without your fingers/latex gloves ever coming into contact with the dressing and wound. This greatly minimizes the risk of infection and avoids the waste of dressings.

Easy removal with minimal skin damage or irritation

Biatain Silicone has a silicone adhesive that combines the softness of a gel with the strength of an adhesive. Silicone is gentle on the skin, making the removal of the dressing far less painful for your patient than traditional adhesives. And it can be reapplied without losing the strength of the adhesive.



Ultra flexible design for superior fit to wound and body

Biatain Silicone has only two layers, making the dressing extremely flexible.


Greater flexibility gives better fit – anywhere on the body, even in awkward places such as ankles, elbows and knees. Greater flexibility also means that there's less need to cut to fit – which saves time and eliminates another source of infection.

Full contact foam for improved healing on a wide range of moist wounds

Biatain Silicone has silicone where you need it – and foam where you need.

The full contact foam can expand to fill the wound bed. Its unique 3D polymer structure ensures high retention even under compression therapy. Together, creates the optimal environment for improved healing of moist wounds.

The foam is free to expand to fill up the wound bed. Unique 3D polymer structure ensures high retention even under compression therapy. This creates the optimal environment for improved healing of moist wounds.



Indications

Biatain Silicone is the versatile solution for moist wounds such as leg ulcers, non-infected diabetic foot ulcers, second degree burns, donor sites, post operative wounds, skin abrasions or pressure ulcers.

The dressing should not be left in place for more than 7 days without clinical inspection by a physician.

Quick guide to Biatain® Silicone assortment



Code	Size	Foam size	Units/Box	Reimbursement Codes
33434	7,5 x 7,5	4 x 4	10	Local subs to provide
33436	12,5 x 12,5	8 x 8	10	Local subs to provide
33437	15 x 15	10 x 10	5	Local subs to provide
33438	17,5 x 17,5	12,5 x 12,5	5	Local subs to provide

* Except pressure ulcers with a high degree of friction and shear force. In these cases, Coloplast recommends the use of Biatain Adhesive dressings.

Ostomy Care
Urology & Continence Care
Wound & Skin Care

 **Coloplast**

Coloplast develops products and services that make life easier for people with very personal and private medical conditions. Working closely with the people who use our products, we create solutions that are sensitive to their special needs. We call this intimate healthcare. Our business includes ostomy care, urology and continence care and wound and skin care. We operate globally and employ more than 7,000 people.

Coloplast A/S
Høtveddam 1
3050 Humlebæk
Denmark

The Coloplast logo is a registered trademark of Coloplast A/S. © 2010-01
All rights reserved Coloplast A/S, 3050 Humlebæk, Denmark.

www.coloplast.com